


## Bisphenols by online SPE and LC-MS/MS


Based on the toxicity of bisphenols to human health and aquatic ecosystem, BPA and its analogues should be strictly controlled in different environmental matrices at low concentrations. The sample preparation step is time-consuming and critical for analytical determination of BPs due to complex analytes pre-concentrating. Online SPE-LC-MS/MS method provides sensitive, effective, rapid, and simple extraction and analysis of bisphenols in water sample matrix.




*Bisphenols on ASTRA® C18 AQ column*

	Precursor (m/z)	Product (m/z)	Collision energy (V)
<b>Bisphenol S</b>	249	108	20
<b>Bisphenol F</b>	199	93	15
<b>Bisphenol AF</b>	335	265	14
<b>Bisphenol A</b>	227	212	11
<b>Bisphenol B</b>	241	212	11


## Bisphenols by online SPE and LC-MS/MS

<b>Online SPE Column</b>	ASTRA® C18 AQ, 10 µm			
<b>Dimensions</b>	30 mm × 2.1 mm			
<b>Part number</b>	AST-5832-PD21			
<b>Loading time</b>	0–3 min			
<b>Analytical Column</b>	ASTRA® C18 AQ, 5 µm			
<b>Dimensions</b>	50 mm × 2.1 mm			
<b>Part number</b>	AST-5832-LG21			
<b>Elution time</b>	3–20 min			
<b>Mobile phase</b>	A: DDW B: MeOH			
<b>Gradient elution</b>	<b>Time (min)</b>	<b>A (%)</b>	<b>B (%)</b>	<b>Flow rate (µL/min)</b>
	0	100	0	450
	1	100	0	1500
	3	100	0	1500
	3.05	100	0	450
	6	100	0	450
	10	0	100	450
	15	0	100	450
	15.05	100	0	450
	20	100	0	450
<b>Temperature</b>	23 °C			
<b>Injection volume</b>	1 mL			
<b>Detection</b>	online SPE-LC-MS/MS-negative MRM mode			
<b>MS instrument</b>	Bruker EVOQ® DART-TQ+ 			
<b>Sample</b>	5 ng/L in DDW			
<b>Analytes</b>	<b>1. Bisphenol S, CAS number 80-09-1</b> <b>2. Bisphenol F, CAS number 620-92-8</b> <b>3. Bisphenol A, CAS number 80-05-7</b> <b>4. Bisphenol B, CAS number 77-40-7</b> <b>5. Bisphenol AF, CAS number 1478-61-1</b>			

## Bisphenols by online SPE and LC-MS/MS




*Bisphenol S*




*Bisphenol F*



*Bisphenol A*



*Bisphenol B*



*Bisphenol AF*